smartLAB[®] move

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smartLAB[®] move

3-D Step-counter with wireless data transfer

Quick Start Instruction

Please read this instruction thoroughly before first using this device


HMM Diagnostics GmbH
D-69221 Dossenheim, Germany
www.hmm.info

Before the first use

1. Remove the battery lash

To turn on the device, you need to remove the lash on the side of the device (pict 1). Please make sure that your gateway is turned on. Please remove the battery lash carefully. After removal the device will turn on. The three LEDs will light up after each other several times for 5 seconds. Then the LEDs will start blinking which indicates that the device is trying to communicate with the gateway to set the device.


pict 1

2. Mounting the carrying cord

You can wear your **smartLAB® move** with the carrying cord included in the set. To do so, you will have to connect the cord to the device by doing the following steps:

Push the thin end of the cord through the two tiny eyes at the back side of the device (pict 2). Pushing together the cord a little bit makes it easier to squeeze it through them (pic. 3). After that, push the thicker end of the cord through the loop on the thin end and pull it tight (pic. 4).

pict 2

pict 3

pict 4

3. Exchanging the Batteries

To exchange the batterie you need to open the battery lid with a coin in direction of the arrow (pict 5). Remove the battery carefully and insert the new battery into the battery lid with the plus side to the lid. Then close the lid by using the coin again.

pict 5

Regularity and Safety Notice

Declaration of Conformity

This product is compliant with the EC Directives: -1999/5/EC

Safety Statement

This product has been verified to conform to the safety certification in accordance with the specification of EN Standards: EN 60950-1:2006 + A11:2009 + A1:2010 + A12:2011

WEEE note

The WEEE (Waste Electrical and Electronic Equipment) Directive, that came into effect as European law on 13th February 2003, led to a major change in the disposal of electrical equipment.

The primarily purpose of this directive is to prevent of electronic waste (WEEE), while encouraging for reuse, recycling and other forms of repro-

cessing to reduce waste.

The WEEE logo on the product and the packaging indicates that it is not allowed to dispose of the product in the ordinary household waste. It is your responsibility to deliver all the disused electrical and electronic equipment to the respective collection points. A separate collection and proper recycling of electronic waste helps dealing economically with natural resources. Furthermore, the recycling of electrical waste is a contribution to the conservation of our environment and thus the health of all people.

More information concerning the disposal of electrical and electronic devices, reprocessing and the collection points

you can get from the local authorities, waste disposal companies, from retailers and manufacturers of the device.

RoHS compliance

This product complies with Directive 2002/95/EC of the European Parliament and the Council of 27th January 2003 regarding the limited use of dangerous substances in electrical and electronic equipment (RoHS) and its variations.

Warranty

HMM Diagnostics GmbH products need to fulfill high quality requirements. Because of this reason, HMM Diagnostics GmbH gives a 2-year warranty by purchasing this **smartLAB**[®] product. You can even extend the warranty from 3 to 5 years without extra pay when you register your product. Please register online under following website: www.hmm.info/en/registration.

Wear parts, batteries etc. are excluded from warranty.

Hersteller/Manufacturer:

HMM Diagnostics GmbH
Friedrichstr. 89
D-69221 Dossenheim, Germany

E-mail: info@hmm.info
www.hmm.info

Weitere Informationen zu den smartLAB® Produkten /
More information on our smartLAB® products:

www.smartlab.org